

Disposable Monopolar EMG Needle Electrode - trocar point



At Technomed we understand the importance of a crystal clear signal, ease of use and patient comfort.

Our new trocar point Monopolar Needle Electrode and cable solutions are designed to give you optimal sharpness and ease of use.

The new hub and cable provide a firm grip and optimal precision in needle positioning.

Optional noise reduction cable available.

User Benefits

- Clear and reliable signal with new recording area design
- Optimal ease of penetration and repositioning with new sharper needle tip
- Precise handling with new ergonomic hub and cable design
- Approved for stimulation
- Optional 50/60 Hz interference elimination with shielded cable

Pre-wired

Needle Length	Diameter	Hub Color	Part Number
25 mm	0.35 mm (28 Gauge)	Red	TE/AP2535-335/01
25 mm	0.45 mm (26 Gauge)	Yellow	TE/AP2545-335/01
37 mm	0.35 mm (28 Gauge)	Orange	TE/AP3735-335/01
37 mm	0.45 mm (26 Gauge)	Green	TE/AP3745-335/01
50 mm	0.45 mm (26 Gauge)	Blue	TE/AP5045-335/01
75 mm	0.45 mm (26 Gauge)	Violet	TE/AP7545-335/01

Detachable

Needle Length	Diameter	Hub Color	Part Number
25 mm	0.35 mm (28 Gauge)	Red	TE/AD2535-008/01
25 mm	0.45 mm (26 Gauge)	Yellow	TE/AD2545-008/01
37 mm	0.35 mm (28 Gauge)	Orange	TE/AD3735-008/01
37 mm	0.45 mm (26 Gauge)	Green	TE/AD3745-008/01
50 mm	0.45 mm (26 Gauge)	Blue	TE/AD5045-008/01
75 mm	0.45 mm (26 Gauge)	Violet	TE/AD7545-008/01

Shelf life for unopened pouches is 36 months
25 individually packed electrodes per box
EtO sterilized



Our Reusable Cables for our Detachable Monopolar Needle Electrodes have a perfect fit to our needles and transmit a reliable and qualitative signal.


Cable Type	Wire Length	Part Number
Shielded	1 m	TE/D54-451
Unshielded	0.6 m	TE/D34-231

Specifications

Properties

Lead-free	✓
Latex-free	✓
BPA-free lead wire	✓
DEHP-free lead wire	✓
RoHS compliant	✓
Lead wire core	Tin plated Copper
Lead wire insulation	PVC
Connector	DIN 42802


Trocar point



USA
Rx ONLY



Your Technomed distributor:



Technomed
medical accessories

Taking you further. Step by step.

Technomed Europe
Amerikalaan 71
6199 AE Maastricht-Airport
The Netherlands
T +31 43 408 68 68
info@technomed.nl
www.technomed.nl