


Disposable Adhesive Surface Electrode



At Technomed we understand the importance of a crystal clear signal, ease of use, repositionability and good grounding.

Our new Disposable Adhesive Surface Electrode is designed with an Ag/AgCl sensor and an extra thick and flexible conductive gel. The electrodes, for recording, stimulation and ground use, are highly flexible and can be repositioned on all types of skin.

The electrode is designed for patient comfort and user friendliness with soft wires and an optimized connector for easy insertion and extraction.

User Benefits

- Clear and reliable signal with optimized sensor of frequency stable Ag/AgCl
- Optimal patient comfort with skin friendly soft and flexible materials
- Extra large soft and flexible gel with high repositionability
- Optimized connector for easy headbox insertion and extraction

Available variants

Disposable 4-Disk Electrodes	Packing	Part Number
20 mm diameter disk x 4; 1 m lead wires: red, blue, green, black	10 sets	TE/K50430-001
20 mm diameter disk x 4; 2 m lead wires: red, blue, green, black	10 sets	TE/K50430-002

Disposable Ground Electrode	Packing	Part Number
40 x 50 mm surface; 1 m lead wire	20 pieces	TE/K50431-002
40 x 50 mm surface; 2 m lead wire	20 pieces	TE/K50431-003

Shelf life for unopened pouches is 24 months

Specifications

Properties	
Lead-free	√
Latex-free	√
BPA-free lead wire	√
DEHP-free lead wire	√
RoHS compliant	√
Lead wire core	Tin plated Copper
Lead wire insulation	PVC
Connector	DIN 42802

USA
Rx ONLY



Your Technomed distributor:




Technomed
 medical accessories

Taking you further. Step by step.

Technomed Europe
 Amerikalaan 71
 6199 AE Maastricht-Airport
 The Netherlands
 T +31 43 408 68 68
 info@technomed.nl
 www.technomed.nl