


# Disposable EEG Cup Electrode Silver/Silver Chloride



**Technomed understands the importance of clear signals, efficiency and ease of use.**

**Technomed's Disposable Ag/AgCl EEG Cup Electrode is designed to optimize EEG procedures to enable focusing on the patient and the EEG recording.**

**The electrode is designed with flexible wires and optimized connector for easy headbox insertion and extraction.**


#### **User Benefits**

- No cleaning and preparation of EEG electrodes needed
- Clear and reliable signal due to optimized layer of frequency stable Ag/AgCl
- Optimized signal quality due to larger electrode surface
- Short cup-to-cable connection to reduce lift risk
- Optimized connector for easy headbox insertion and extraction
- Optimal patient comfort with low profile disk, soft and flexible anti-tangle lead wires
- Reduced risk of cross-contamination

## Available variants

Materials	Cup Size	Cup Height	Wire Length	Packing	Part Number
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1 m	Sets of 25, 5 x 5 colors	TE/C52-438
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1.5 m	Sets of 25, 5 x 5 colors	TE/C52-638
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	2 m	Sets of 25, 5 x 5 colors	TE/C52-838
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	2.5 m	Sets of 25, 5 x 5 colors	TE/C52-938
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1m	Sets of 10, 2 x 5 colors	TE/C52-432
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	1.5 m	Sets of 10, 2 x 5 colors	TE/C52-632
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	2 m	Sets of 10, 2 x 5 colors	TE/C52-832
Silver/Silver Chloride (Ag/AgCl)	10 mm	2.6 mm	2.5 m	Sets of 10, 2 x 5 colors	TE/C52-932

**5 Colors:** red, green, white, yellow, blue



## Specifications

Cup	AgCl layer
Lead-free	√
Latex-free	√
BPA-free	√
DEHP < 0.1% by weight	√
REACH & RoHS compliant	√
Lead wire core	Tin Plated Copper
Lead wire insulation	PVC
Connector	DIN 42802

Your Technomed Business Partner:

USA
Rx ONLY



  
**Technomed**  
 medical accessories

Taking you further. Step by step.

### Technomed

Amerikalaan 71  
6199 AE | Maastricht-Airport | The Netherlands  
T +31 43 408 68 68 | info@technomed.nl

### Technomed USA

Suite 4990 Constellation Drive  
White Bear Township | MN 55127 | United States  
T +1 651 288 0890 | info@technomed-usa.com

[www.technomed.nl](http://www.technomed.nl)