

# Reusable EEG Cup Electrode Gold and Silver



**Technomed understands the importance of a crystal clear signal, ease of use and electrode durability.**

**Technomed's Reusable EEG Cup Electrode Gold and Silver are designed to provide reliable signals and longtime ease of use.**

**These electrodes are designed for easy adjustment to patient contours to avoid electrode lift.**

**The connector is optimized for easy headbox insertion and extraction.**

#### **User Benefits**

- Clear and reliable signal.
- Short cup-to-cable connection to reduce lift risk.
- Optimal patient comfort with low profile disk, soft and flexible anti-tangle lead wires.
- Optimized connector design for easy headbox insertion and extraction.
- Long electrode life with watertight design and Kevlar reinforced wires.
- Available in two different cup sizes.

## Available variants

Materials	Cup Size	Cup Height	Wire Length	Packing	Part Number
Silver Plated (Ag)	10 mm	3.1mm	1.0 m	Sets of 12, 2 x 6 colors	TE/C12-434
Silver Plated (Ag)	10 mm	3.1 mm	1.5 m	Sets of 12, 2 x 6 colors	TE/C12-634
Silver Plated (Ag)	10 mm	3.1 mm	2.0 m	Sets of 12, 2 x 6 colors	TE/C12-834
Silver Plated (Ag)	10 mm	3.1 mm	2.5 m	Sets of 12, 2 x 6 colors	TE/C12-934
Silver Plated (Ag)	6 mm	1.7 mm	1.0 m	Sets of 12, 2 x 6 colors	TE/C14-434
Silver Plated (Ag)	6 mm	1.7 mm	1.5 m	Sets of 12, 2 x 6 colors	TE/C14-634
Silver Plated (Ag)	6 mm	1.7 mm	2.0 m	Sets of 12, 2 x 6 colors	TE/C14-834
Gold Plated (Au)	10 mm	3.1 mm	1.0 m	Sets of 12, 2 x 6 colors	TE/C32-434
Gold Plated (Au)	10 mm	3.1 mm	1.5 m	Sets of 12, 2 x 6 colors	TE/C32-634
Gold Plated (Au)	10 mm	3.1 mm	2.0 m	Sets of 12, 2 x 6 colors	TE/C32-834
Gold Plated (Au)	10 mm	3.1 mm	2.5 m	Sets of 12, 2 x 6 colors	TE/C32-934
Gold Plated (Au)	6 mm	1.7 mm	1.0 m	Sets of 12, 2 x 6 colors	TE/C34-434
Gold Plated (Au)	6 mm	1.7 mm	1.5 m	Sets of 12, 2 x 6 colors	TE/C34-634
Gold Plated (Au)	6 mm	1.7 mm	2.0 m	Sets of 12, 2 x 6 colors	TE/C34-834

**2 x 6 colors:** red, green, white, yellow, blue, brown



## Specifications

Cup	Ag or Au plating
Not made with natural rubber latex	√
REACH & RoHS (EU) compliant	√
Lead wire design	Tinsel wire
Lead wire core	Kevlar and Silver Plated Copper
Lead wire insulation	TPU
Spring loaded connector	Gold Plated Brass / PVC
Connector	DIN 42802-1

Your Technomed Business Partner:



**Cleaning and disinfection:** see IFU for details



**Legal manufacturer:**

**Technomed Europe**

Amerikalaan 71

6199 AE | Maastricht-Airport | The Netherlands

T +31 43 408 68 68 | info@technomed.nl

**Technomed USA**

17200 Medina Road, Suite 800

Plymouth | MN 55447 | United States

T +1 651 288 0890 | info@technomed-usa.com

[www.technomed.nl](http://www.technomed.nl)