

Disposable Stimulation Clip



Technomed understands the importance of clear signals, efficiency and ease of use consumables.

The Disposable Stimulation Clip turns a metal surgical dilator into a monopolar stimulating probe.

Fits and connects to most common dilators used in spine surgery.

High quality stimulus conduction when applied to surgical dilators.

User Benefits

- Suitable for EMG, NCS, IONM and EP monitoring.
- Easy to attach to dilators, with a diameter between 5mm and 22mm.
- Insulated ergonomic handle.
- Standard 1.5mm DIN 42802-1 connector with a 2.5m lead wire that can connect to a wide range of nerve monitors.
- Ready to use, no cleaning or preparation needed.
- Individually sterile packed.
- Single patient use to avoid cross-contamination.

Disposable Stimulation Clip

Materials	Width of Clip	PVC Mold	Wire Length	Packing	Part Number
Stainless Steel	6mm	Black	2.5 m	Box of 1	TE/Q95-931



Specifications

Not made with natural rubber latex	√
REACH & RoHS (EU) compliant	√
Lead wire core	Tin Plated Copper
Lead wire insulation	PVC
Connector	DIN 42802-1
Shelf life for unopened pouches	36 months

Your Technomed Business Partner:

USA
Rx ONLY

CE
0344


Technomed
medical accessories
EST. 1980


40
YEARS
OF MEDICAL
INNOVATIONS

Legal manufacturer:

Technomed Europe

Amerikalaan 71

6199 AE | Maastricht-Airport | The Netherlands

T +31 43 408 68 68 | info@technomed.nl

Technomed USA

17200 Medina Road, Suite 800

Plymouth | MN 55447 | United States

T +1 651 288 0890 | info@technomed-usa.com

www.technomed.nl